

D:O IMPLANT

Dental X-Ray Machine

HIGH LENTIFON DESIGN

Optimalisasi desain
konstruksi untuk meningkatkan
efisiensi radiasi dan mengurangi
paparan radiasi.

Desain ergonomis yang memudahkan
penggunaan dan mengurangi beban
fisik operator. Fitur keamanan yang
meningkatkan keselamatan pasien dan
operator.

Info & Kontak

021-2251-1234



Dental X-Ray Machine



The Handheld Dental X-Ray machine is a cutting-edge choice for dental diagnostics, featuring ergonomic design, long battery life, and intelligent angle adjustment for precise imaging.



Long battery life



Smart camera angle assistance feature



Ergonomic design

PARAMETERS INSTRUCTION

	Model	POCT22		Anode Angle	12.5 °
	Tube Voltage	70KV±10%		Focus	0.4mm
	Tube Current	2mA±20%		Exposure Time Range	0.04s -2.0s
	Input Voltage	100V~240V		Leak Emissivity	1 meter 0.25mGy/h
	Frequency	50Hz-60Hz		Machine Weight	1.7kg
	X-Tube	Canon (Toshiba)		Relative Humidity	30%RH-75%RH



Compact & Portable



User Friendly Interface



High Resolution Imaging



Stable & Reliable



D:O IMPLANT

DIO Digital Implant India Pvt Ltd
706-711, 7th Floor, Vipul Square,
Sushant Lok Phase-1, Sector-28,
Gurugram-122009

Tel. +91 0124-4751803 | 7303331511
dioimplant.co.in

